

Velino Care

08557 The Velino nursing bed is a multifunctional height-adjustable column bed that maximizes both the patient comfort and the user-friendliness.

Carriage

- High/low adjustment by means of 2 telescopic lifting columns.
- Reverse Trendelenburg adjustment possible up to $\pm 14^\circ$.
- The castors are hidden by means of a synthetic cover cap with a wooden look or with a grey uni colour.
- Twin castors \varnothing 50 mm, braked 2 per 2.
- The design brake pedal can be activated on both sides of the bed.
- The open carriage is very accessible and easy to clean.

Mattress platform

- Three- or four- (*) section mattress platform.
- The backrest is electrically adjustable up to $\pm 70^\circ$ with autoregression.
- The knee rest is electrically adjustable up to $\pm 34^\circ$ (in case of 4 sections).
- The footrest is manually adjustable for the Fowler position (in case of 4 sections).
- The legrest is manually adjustable up to $\pm 14^\circ$ (in case of 3 sections).
- The electric mattress platform manipulations are equipped with a clamp safety.
- The mattress platform is made of rounded steel profiles.
- The mattress platform sections are equipped with openings for patient restraint straps.
- The backrest and footrest have integrated mattress holders.
- Telescopic bed extension possible up to 20 cm (not with Pollux sideguards).
- The head end is equipped with 2 cases to install an IV rod or a lifting pole.
- No visible fastening bolts in the foot panel.
- The bed panels are made of solid oven-dried beech (*).
- Rounded corners to increase the passive safety.

Electric bed manipulations

- All electric bed manipulations are activated by means of maintenance-free Linak actuators.
- The mattress platform actuators are standard equipped with a clamp safety and can be replaced without needing any tools.
- The electric bed functions can be activated by means of a hand switch or an *mini-nurse box* (*).
- The hand switch has a preprogrammed sleeping and relax function.
- A *mini-nurse box* (*) makes it possible to switch off several functions selectively, whereas a lock box (*) switches these functions off by means of 1 single button.

Paneltypes

- Midea (standard)
- Valentino (*)
- Puro (*)
- Felina (*)

Sideguards

- Pollux (*)
- Trix (*)


Accessories

We would like to refer to our accessory list.



Dimensions and technical data

- Width: 100,5 cm. Length: 208 cm
- Height variation of the mattress platform: $\pm 40,5$ cm - $\pm 80,5$ cm
- Safe working load: 250 kg
- Weight (without options): 140 kg
- Mattress platform dimensions: 204 x 86 cm
- Mattress dimensions: 195 x 85 cm Thickness: 12 to 15 cm
- Standby consumption: 0,5 W (max. 4,38 kWh/year)



Materials and finishing

- Wood: beech (massive, veneer), chipboard (melaminated).
- Metal: epoxy-coating RAL 9006 white aluminium.
- Self-lubricating synthetic hinges (= maintenance-free).
- Withstands the common cleaning products.

For more information about the **materials** and **construction**: please consult our material lists. We would refer to our extensive **colour gamut**.

Options

- 4-section mattress platform.
- 3- or 4- section mattress platform with removable solid laminate sections.
- Foot extension up to + 20 cm.
- Panel types: Puro, Felina or Valentino.
- Bumpers (not possible on Midea).
- Wall spacer.
- Sideguards.
- Total lock device (ACL).
- Mini-nurse box (ACC).
- Hand switch with selective locking system.
- Nightlight under the bed.
- Accessory hooks.
- Wooden surrounds of the mattress platform frame.
- Battery.
- Manual quick release of the backrest (CPR).
- Extendible drawer

(*) Option



HAELVOET NV | L. BEKAERTSTRAAT 8, 8770 INGELMUNSTER, BELGIUM | T +32 (0) 51 48 66 95 | INFO@HAELVOET.BE | WWW.HAELVOET.BE

Since Haelvoet continuously aims to make technological progress, Haelvoet reserves the right to modify the product.